

G.E.M.M.® in Practice: Clinician Mentoring Programme

ACCREDITATION PATHWAY
CLINICIAN MENTORING
PRACTICE-READY



From Foundations to Accreditation - Practice-ready implementation

The G.E.M.M.® Mentoring Programme supports practitioners as they implement the G.E.M.M.® Protocol with confidence and clinical precision, bridging the gap between theory and application.

This structured, three-part program guides you step-by-step through the Foundations Course and each upskill session, through to the final Accreditation Assessment, ensuring understanding of core concepts and best-practice clinical implementation.

Throughout the programme, you will develop:

- Stronger understanding of core nutrigenomic principles in practice
- Greater confidence in formulation selection
- Safe and effective titration of sulforaphane
- Practical implementation strategies across common patient types
- Applied clinical thinking through live case discussion



[LEARN MORE ABOUT G.E.M.M.®](#)

WAS \$495 INC GST

NOW \$395 INC GST

PLACES CAPPED AT 50 PRACTITIONERS

WHAT'S INCLUDED:

- 18.75 hours CPE across Foundations & Upskill sessions
- 3 live mentoring sessions with the Cell-Logic Education Team
- Level A Starter Pack:
GFD 120c, MCP-U 250g, Gemmune I.B 56c valued at \$118 inc GST wholesale price
- Addition to Cell-Logic's preferred clinician online directory
- G.E.M.M.® marketing kit for digital platforms
- Support towards G.E.M.M.® Accreditation

June 2026 Live Mentoring Series

12pm to 1pm AWST

Thursday 11 June 2026 Core principles, implementation foundations, product advisory, cautions and contraindications


Thursday 25 June 2026 Best-practice implementation, patient types, product selection, tracking and testing

Thursday 9 July 2026 Case submission review, peer discussion and final mentoring

Ready to join us in June?

1. Scan QR code to apply or login
2. Add the G.E.M.M.® Foundations Course to cart.
3. Apply code GEMM-RENER at checkout.

Offers ends 11th June 2026 | Places capped at 50



Phone (+61 7) 3488 0385 or visit cell-logic.com.au